

# Contents

CNV 分析参数	1
CNV 解析度	2

## SNV+CNV

```
dragen -f -r ${1} \  
-1 ${2} -2 ${3} --output-directory ${4} --output-file-prefix ${5} \  
--RGID Illumina_RGID --RGSM ${5} \  
--enable-map-align true \  
--enable-map-align-output true \  
--output-format BAM \  
--enable-duplicate-marking true \  
--enable-sort true \  
--enable-variant-caller true \  
--vc-enable-vcf-output true \  
--enable-vcf-compression true \  
--enable-cnv true \  
--cnv-enable-self-normalization true \  

```

## SV

```
--enable-sv true
```

## repeat

```
--repeat-genotype-enable true  
--repeat-genotype-specs /opt/edico/repeat-specs/hg19_expanded/
```

## gvcf

```
--vc-emit-ref-confidence GVCF
```

**DRAGEN-ML:** typically removes 30-50% of SNP FPs, with smaller gains on INDELS. FN counts are reduced by 10% or more.

```
--vc-ml-dir=/opt/edico/resources/ml_model/hg19 --vc-ml-enable-recalibration=true
```

**Star Allele Caller** 星号等位基因：用于药物基因组

```
--enable-star-allele true
```

**Targeted callers:** 输出.targeted.vcf.gz 和 .targeted.json

```
--enable-targeted=true
```

## High Sensitivity Mode

```
--vc-enable-high-sensitivity-mode=true
```

## Pharmacogenomics (PGx) grouping

```
--enable-pgx=true
```

# CNV 分析参数

```
--enable-cnv true  
--cnv-filter-copy-ratio 0.2 # The default value is 0.2, leading to calls less than CR=0.8 or  
↪ greater than CR=1.2.  
--cnv-filter-length 10000 # Specifies the minimum event length in bases at which a reported  
↪ event is marked as PASS in the output VCF file. The default is 10000  
--cnv-filter-qual 10 # PASS in the output VCF file
```

## CNV 解析度

WGS_Coverag_per_Sample	Recommended_Resolution(bp)
5X	10000
10X	5000
>=30X	1000

*-cnv-interval-width* 用来控制解析度，WES 默认是 500，WGS 默认是 1000 该参数在分析是需要设置，如果设置变小会增加分析时间